SINCE 1966
MADE IN ITALY

Smarty

by **FLAEM**®

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

ACO567P

Smarty

**Portable device
for aerosol therapy using
the latest in PMVT technology (Passive MESH Vibrating Technology)**


Thank you for your purchase. Our goal is the complete satisfaction of our customers by offering cutting-edge products for the prevention and treatment of the respiratory tract ailments.


Visit our website www.flaem.it to view the whole range of Flaem's products.

**Read these instructions carefully
in order to use the device correctly.
We recommend that you store this manual
for future consultation.
Smarty is a medical device for home use.**

Intended use: this is an electronic Medical Device for aerosol therapy using PMTV technology (Passive MESH Vibrating Technology). **Smarty** is suitable for nebulising **medications which are either in solution** eg mucolytics, or **suspension** eg corticosteroids such as **beclomethasone dipropionate** and **budesonide**) and saline solution.

WARNING
TO PREVENT DAMAGE TO THE MESH,
CAREFULLY FOLLOW STEPS 1 TO 5 AS LISTED BELOW

LIST OF DEVICE COMPONENTS

1 - Main Unit


- 1A - ON DEMAND nebulisation button
- 1B - CONTINUOUS nebulisation / START AND STOP button
- 1C - Medication chamber release button
- 1D - Battery Compartment Cover
- 1E - Battery Compartment

2 - Removable medication chamber

- 2A - Medicine cup
- 2B - Nebulising chamber
- 2C - Medication cup locking clip
- 2D - Medication cup cover
- 2E - Nebulising chamber and Medication cup release clips
- 2F - Accessories connection port
- 2G - Nebuliser head

Accessories

- 3 - Mouthpiece
- 4 - Child mask
- 5 - Adult mask
- 6 - 2 1.5V AA alkaline batteries
- 7 - Carrying bag

IMPORTANT WARNINGS


- As with any electrical device, especially in the presence of children, the use of **Smarty** requires some basic safety precautions.
- The device should not be used as a toy. Please take particular care when it is used by children.
- If the device is used by or in the presence of children or people who require special assistance, the supervision of an adult who has read this manual is necessary.
- Some parts of the device are small enough to be swallowed by children; keep the device out of reach of children.
- Never use batteries other than those recommended in this manual.
- Use only batteries of known brands. Do not use old and new batteries.
- AA alkaline or Lithium batteries rechargeable Ni-Mh can be used. Rechargeable batteries shall be charged using a battery charging device that is not supplied with this device.
- The average battery life depends on the brand used.
- It is very important to use batteries of the same brand and same type.
- Remove the batteries before storing the unit when it is not used for long periods.
- Smarty should not be used in the presence of an anaesthetic mixture which is flammable with oxygen or nitrous oxide.
- **The device housing has been designed for limited contact with liquids (no immersion). Do not wash the device under running water or by immersion. Carefully follow the "CLEANING, SANITISATION, DISINFECTION" in this manual. Absolutely do not use alcohol to clean the device**
- Do not use the device while taking a shower or a bath.
- Do not operate the device while driving or in any other situation where a distraction could be dangerous for the user, for close people or animals or surrounding objects.
- Do not expose the device and the batteries to extreme temperatures. Keep the device and batteries away from heat sources, direct sunlight or excessively hot and humid environment.
- The average expected life of the device is approximately 300 hours.
- Repairs must be performed by authorised FLAEM personnel only. Any unauthorized repairs will void the warranty and may pose a safety hazard for the user.
- **WARNING!** Do not modify this device without authorisation of the manufacturer.
- The Manufacturer, Importer and Seller will be held responsible for the safety, reliability and performance, only if the device is used in accordance with the instruction manual.
- **Microbial contamination:** personal use of the device and accessories is recommended to avoid any risk of infection.
- **Interactions:** the materials used have been tested according to the standards

of biocompatibility (ISO 10993-5 and ISO 10993-10) in compliance with the essential requirements of the 93/42 EEC Medical Device Directive.

- If the device is used with different medications, it is recommended to remove all residues. Therefore, clean it after each treatment so as to obtain the highest degree of hygiene and to optimise the life and operation of the device
- The device does not contain parts to be repaired by the user. The warranty does not cover batteries or damage caused by incorrect, used up or improperly stored batteries.
- The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation.

OPERATING INSTRUCTIONS


1 - BATTERY INSTALLATION

1.a Use a small coin to open the **Battery Compartment cover (1D)** located on the bottom of the unit.

1.b Insert the two **AA batteries (6)**, in the **Battery Compartment (1E)**, making sure that the polarity is correct.

1.c Close the **Battery Compartment Cover (1D)** making sure that the cover is properly.


2 - PREPARATION

WARNING!

Before each use, the unit and the accessories should be cleaned and/or disinfected according to the **CLEANING, SANITISATION, DISINFECTION** instructions (paragraph 6)

We recommend personal use, both of the nebuliser and accessories to avoid possible risk of infection. The average life span of the accessories is approximately 1 year.

Use only FLAEM original accessories.

Wash your hands thoroughly before preparing the medication

2.1 Lift the **locking clip (2C)** of the **Medication cup cover (2D)**

2.2 Open the **Medication cup cover (2D)** completely


2.3 Pour the medicine (MAX 8 ml) into the **Medication cup (2A)**

2.4 Close the **Medication cup cover (2D)**

2.5 Fasten the **locking clip (2C)**

WARNING!

After pouring the medicine into the **Medication cup (2A)**, if the **cover (2D)** is open **DO NOT** press the **release clips (2E)**, as the **Nebulising chamber (2B)** might open causing the medicine to spill.


WARNING!

The medicine must be nebulised within 10 minutes of being put in the **Medication cup (2A)**

3 - UNIT CONFIGURATION




3.1 Connect the desired accessory to the **Accessories connection port (2F)**:

- **Mouthpiece (3)**
- **Child mask (4)**
- **Adult mask (5)**

If a mask is used, place it on the face as shown in the picture

4 - OPERATION

Select the preferred delivery mode ("**ON DEMAND**" or "**CONTINUOUS**")

4.1 ON DEMAND button (1A) Press this button when inhaling to obtain the delivery of medication "**ON DEMAND**". The device works only when you press the button.

4.2 START AND STOP / CONTINUOUS button (1B): Press this button to obtain the "**CONTINUOUS**" delivery of medication (fast treatment). To switch off the unit, press the **button (1B)** again.

4.3 When the unit is operating "**ON DEMAND**", the **button (1A)** lights up in **BLUE** colour.

4.4 During operation in "**CONTINUOUS**" mode, the button (1B) lights up as follows:

4.4.a STEADY BLUE: the device is operating normally

4.4.b FLASHING RED: batteries are used up; replace with new batteries of the same type

4.5 Sit comfortably holding the device in your hand, place the **Mouthpiece (3)** in your mouth then breath in and out slowly through the mouth. Alternatively, you can use the **Child Mask (4)** or the **Adult Mask (5)**.

4.6 To increase the effectiveness of treatment, after inhaling hold your breath for a moment, so that the inhaled medication can deposit along the respiratory tract. Then exhale slowly.

4.7 The unit turns off automatically:


4.7.a after 5 minutes of use

4.7.b when, during use, the **medication chamber (2)** is accidentally removed .


5 - DISASSEMBLY

When treatment is finished:



- 5.1 Press the **Medication chamber release button (1C)** positioned on the back of the unit.
- 5.2 Remove the **Medication chamber (2)**.
- 5.3 Carefully clean it according to the **CLEANING, SANITISATION, DISINFECTION** instructions in paragraph 6.

- 5.4 To access the **Nebuliser head (2G)** in the **Nebulising chamber (2B)** release the **Medication cup locking clip** and open the **Medication cup cover (2D)**.
- 5.5 Press simultaneously the two **release clips (2E)** on both sides of the Medication chamber.
- 5.6 Rotate until the **Nebulising chamber (2B)** opens completely.


- 5.7 Use your finger to pull out the **Nebuliser head (2G)** from the bottom of the chamber.

WARNING!


- Do not use tools: this may cause irreparable damage to the **Nebuliser head (2G)**.
- Do not press the centre of the **Nebuliser head (2G)** with your fingers or tools from the inside: this may cause irreparable damage.

6 - CLEANING, SANITISATION, DISINFECTION

Before and after each use, the device and the accessories must be cleaned and disinfected properly as described below. If this is not done, some microorganisms may collect in the device, causing the risk of infection. Absolutely do not use alcohol to clean the device

6.1. CLEANING Accessories and Removable Upper Part (2):


6.1.1. WASHING: Wash the components **2-2G-3-4-5** with hot drinking water with mild detergent for dishes (non-abrasive), then rinse thoroughly with a jet of hot drinking water to remove any detergent residue.

6.1.2. DRYING

6.1.2.a. After washing and rinsing all the components, shake them to remove excess water and reassemble them for the time being, as described in the following paragraph **7. REASSEMBLY OF THE DEVICE**, but without the application of the **Mouthpiece (3)** or **Masks (4)** and **(5)**. At this point, to eliminate the residual water that has deposited in the **nebuliser head (2G)**, press **Button 1A** for about 2-3 seconds until the device stops spraying. This operation is very important, because it prevents the deposit of limestone in the micropores of the **Nebuliser Head (2G)** which could affect proper nebulisation of the drug in subsequent therapeutic applications.

6.1.2.b. To complete the drying of all the components, disassemble them again and place them on a paper towel as shown in the image. Leave them out and in a dry place (e.g. not in the bathroom) for about 3-4 hours.


ATTENTION:

Rest the **Head of the nebuliser (2G)** keeping the side with the coloured seal turned upwards.

WARNING!

Do not clean the **Mesh (2G)** using brushes or other instruments: this may cause irreparable damage.

6.2 CLEANING the Main Unit


To clean the **Main Unit (1)** use a soft sponge or a cloth slightly moistened with an antibacterial detergent, then dry it with a soft paper towel.

ATTENZIONE!


Do not immerse the **Main Unit (1)** in water and do not wash it under running water.

6.3. SANITISATION

6.3.1. Immerse the opened **Removable Medication Chamber (2)** and the accessories **2G-3-4-5**, in a solution of 50% water and 50% white vinegar, then rinse thoroughly with drinking water.


6.3.2. DRYING: carefully perform all the operations described in the previous paragraph **6.1.2**.



6.4. DISINFECTION

WARNING! In this case, the **Device Body (1)** is not treated. To clean it, please refer to the previous point

6.2. CLEANING Device Body. The parts that can be disinfected are those identified **2G-3-4-5**. The disinfection procedure described in this paragraph is to be performed before and after use of the device and is effective for accessories subjected to this treatment only if all the steps are performed and only if the accessories to be treated are sanitised beforehand as described in the preceding paragraph. The disinfectant must be an electrolytic chloroxidizer (active principle: sodium hypochlorite) specific for disinfecting, which is available in any pharmacy.

6.4.1. Use a container of a suitable size to contain all the components to be disinfected. Fill it with a solution of potable water and electrolytic chloroxidizer (active ingredient sodium hypochlorite) suitable for

disinfection according to the dilution proportions recommended in the disinfectant instructions.


6.4.2. Submerge every component completely in the solution; try to prevent that air bubbles form in contact with the components. Keep the components immersed for the period recommended on the disinfectant instructions based on the selected concentration / dilution.

6.4.3. Remove the disinfected components and rinse them thoroughly under warm potable water.


6.4.4. DRYING: carefully perform all the operations described in the previous paragraph 6.1.2.

6.4.5. Discard solution according to the instructions provided by the disinfectant manufacturer.

7 - DEVICE RE-ASSEMBLY

7.1. Insert the **Mesh(2G)** in its seat located inside the **Front Part (2B)** of the **Removable medication chamber (2)** by grasping the external part with the thumbscrew with your fingers and carefully keeping the part with the coloured seal turned upwards.


7.2. Push the **Mesh (2G)** right down, taking care to always place your fingers on its milled outer part.

WARNING!


Do not press the **Nebuliser head (2G)** in the centre and do not use tools: this may cause irreparable damage.


7.3. Close the **Removable Medication chamber (2)** by rotating the **Medication Cup (2A)** towards the **Nebulising chamber (2B)** until the two **clips (2E)** are locked.

7.4. Close the **Medication Cup Cover (2D)** by means of the **locking clip (2C)**

7.5. Insert the **Removable Medication chamber (2)** onto the **Main Unit (1)** matching the guides on the two parts.

7.6. Push the **Removable Medication chamber (2)** right down until the **Medication container release button (1C)** located in the rear of the **Main Unit (1)** springs with a **CLICK**.

TROUBLE-SHOOTING

PROBLEM	CAUSE	SOLUTION
The device switches off after a few seconds (blue light).	No electrical contact between the Removable Medication chamber (2) and the Main Unit (1)	Check that the Removable Medication chamber (2) is properly hooked to Main Unit (1)
The device switches off	It's been 5 minutes since device switching on	Turn on the unit again
The START AND STOP Button (1B) shows a flashing red light and the device does not nebulise or switches off.	Batteries are used up or low	- Replace the batteries if they are alkaline. - Recharge the Batteries if they are rechargeable Ni-MH or Lithium type.
The device does not switch on.	No electric power supply	- Verify that Batteries are inserted with the correct polarity. - Used up or low Batteries. Replace or recharge.
The device switches on but does not nebulise or nebulises little.	- Malfunction of the Mesh (2G) - No medicine supply in the nebulising area	- Verify that the Nebuliser head (2G) is correctly inserted in its seat and that the Removable Medication Chamber (2) is correctly assembled to the Main Unit (1) as described in paragraph 7 DEVICE RE-ASSEMBLY. - Sanitise the Nebuliser head (2G) as described in paragraph 6. CLEANING, SANITISATION, DISINFECTION. - Verify that the small hole in the centre of the Medication cup cover (2D) is not obstructed. If so clean it with a needle

If after checking the aforementioned conditions, the unit still does not nebulise, we recommend that you contact the dealer or an authorised Flaem service centre.

8 - TECHNICAL CHARACTERISTICS

Smarty Electronic Nebuliser Model P0309EM

Power Supply	3.0Vdc with 2 Batteries AA Alkaline, 1.5Vdc 2.4Vdc with 2 Batteries AA Rechargeable Ni-MH or Lithium, 1.2Vdc
Power	1.5W
Transducer frequency	150 KHz
Noise level (at 1 m)	0 dBA
Device dimensions	51(W) X 67(D) X 140(H) mm
Device weight (excluding batteries)	125g
Carry bag dimensions	140(W) X 80(D) X 150(H) mm
Delivery rate ml/min⁽¹⁾	0.30 (approx)
MMAD⁽²⁾	3.17 µm
Breathable fraction <5µm⁽²⁾	75 %
Medication cup capacity⁽¹⁾	8 ml

⁽¹⁾ Data shown is in accordance with Flaem Nuova internal procedure I29-P07.5 Values shown for the aerosol delivery may vary depending on the patients respiratory Capacity

⁽²⁾ In vitro testing in compliance with the European Standard EN 13544-1 for nebulizing systems. Further details available on request.

APPLIED PARTS

The BF type applied parts are: patient accessories (3,4,5)

9 - OPERATING CONDITIONS

Temperature	min 10°C ~ max 40° C
Air humidity	min 10% ~ max 75% RH
Atmospheric Pressure	min. 690 hPa ~ max. 1060 hPa

10 - STORAGE CONDITIONS

Temperature	min 10°C ~ max 35°C
Air humidity	min 10% ~ max 75% RH
Atmospheric Pressure	min. 690 hPa ~ max. 1060 hPa

11 - SYMBOLS

	Batteries			Battery disposal: Used up batteries must be disposed of in appropriate containers.
	Type BF applied parts			Manufacturer
	WARNING! Check the instructions for use manual			Keep dry
	Direct current			Serial number of the device
	CE Marking medical ref. Dir 93/42 EEC Directive and subsequent updates			

12 - DISPOSAL OF DEVICE

In conformity with Directive 2002/96/EC, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to "sorted waste collection". The user must therefore take (or have taken) the above waste to a pre-sorted waste collection centre set up by the local authorities, or else give it back to the dealer when purchasing a new appliance of the same type. Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects of any incorrect waste management on the environment and public health. The unlawful disposal of the product by the user could result in administrative fines as provided by the laws transposing Directive 2002/96/EC of the European member state or of the country in which the product is disposed of.

13 - ELECTROMAGNETIC COMPATIBILITY


This device has been designed to satisfy requirements currently required for electromagnetic compatibility (EN 60 601-1-2:2007). Electrical medical devices require special care during installation and use with respect to EMC requirements. It is therefore required that they be installed and/or used according to the manufacturer's specification. Potential risk of electromagnetic interference with other devices, in particular with other devices for analysis and treatment. Radio and mobile telecommunications devices or portable RF (mobile phones or wireless connections) may interfere with the operation of electrical medical devices. For further information visit www.flaem.it website. Flaem reserves the right to make technical and functional modifications to the product without notice.

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Expected duration: The duration times of the various components listed below are related to an use of the device which provides for the nebulisation of saline to 2 times a day, 5 minutes each time, at an ambient temperature of 23 °C.
The duration time may vary depending on the environment of use.

Main unit	5 years
Mouthpiece	1 year
Child mask	1 year
Adult mask	1 year
Nebuliser head	6 months
Medication chamber	1 year

AVAILABLE SPARE PARTS AND ACCESSORIES

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
Rivenditore (timbro e firma)

The warranty terms here specified are valid only in Italy for Italian residents. In all other countries, the warranty will be provided by the local dealer that sold you the unit, in accordance with the applicable laws.

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© 2013 FLAEM NUOVA® All rights reserved
FLAEM NUOVA S.p.A.